

N₂ MICRO

	Modello Model	Dimensioni Dimension					Peso Weight	
	Purezza gas % Gas purity %	95%	98%	99%	99.5%	99.9%	(cm.) (Kg.)	
CABINET	MICRO 0.5	1.25	0.80	0.60	0.50	0.30	50x24x75 35	
	MICRO 1	2.50	1.60	1.20	1.00	0.60	45x67x90 55	
	MICRO 2	5.00	3.20	2.40	2.00	1.20		65
	MICRO 3	7.50	4.80	3.60	3.00	1.80		75
	MICRO 4	10.00	6.40	4.80	4.00	2.40		85
	MICRO 5	12.50	8.00	6.00	5.00	3.00	60x60x110 205	
	MICRO 6	16.20	10.40	7.80	6.50	3.90	70x78x185 265	
	MICRO 8	21.20	13.60	10.20	8.50	5.10		315
	MICRO 12	32.50	20.80	15.60	13.00	7.80		400
	MICRO 16	42.50	27.20	20.40	17.00	10.20		450
Coeff. Aria Richiesta Air Demand Coeff.		2.60	3.40	4.00	5.00	7.00		

N.B. - i valori riportati in tabella sono riferiti ad una pressione di ingresso aria di 8 bar ed una pressione di uscita azoto di 6 bar.

- values shown in the table are referred to a pressure of air entry of 8 bar and a pressure of nitrogen exit of 6 bar.